

## Color Options

**Translucent Material - M2000 / M35 / M45 / SlikFit™**  
 • Platinum Cure Silicone • Hypoallergenic • Antibacterial • Comfortable




**Opaque Material - Passguard™ / Softguard™ / Neon Lites™**  
 • Platinum Cure Silicone • Hypoallergenic • Antibacterial • Comfortable



## Color Style Options

### Silicone Earmolds

| COLOR STYLE                                                                                         | DESCRIPTION                                                                             | # COLORS                      | MATERIAL                                                 | GLITTER      | FLUORESCENT  |
|-----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|-------------------------------|----------------------------------------------------------|--------------|--------------|
| <br>Gel-e-Burst  | Clear with up to 4 additional colors swirled in                                         | Clear + up to 4 colors        | M35, M45<br>M2000<br>SlikFit™                            | Added Option | Added Option |
| <br>Rainbow      | 2 to 5 colors, overlapped in a rainbow pattern                                          | 2 to 5 colors                 | M35, M45<br>M2000<br>SlikFit™<br>Passguard<br>Neon-Lite™ | Added Option | Added Option |
| <br>Marble     | White with 1 additional color, forming a marble pattern                                 | White + 1 color               | M35, M45<br>M2000<br>SlikFit™<br>Passguard<br>Neon-Lite™ | Added Option | Added Option |
| <br>Polka-dots | 1 main color with up to 4 additional colors<br><i>1-2 additional colors recommended</i> | 1 main color + up to 4 colors | M35, M45<br>M2000<br>SlikFit™<br>Passguard<br>Neon-Lite™ | Added Option | Added Option |
| <br>Solid      | 1 color of choice                                                                       | 1 color                       | M35, M45<br>M2000<br>SlikFit™<br>Passguard<br>Neon-Lite™ | Added Option | Added Option |