

# LifePlugs®

Take the edge off  
of everyday life.

- MICROSONIC is proud to introduce our new product line, LifePlugs® and LifePlugs® Sleep.

Our new custom fit extended-wear earplugs that tone down the auditory disturbances of life.

**LifePlugs®** Vented & pressure free. Perfect for yoga, exercise, or meditation.

**LifePlugs® Sleep** for a quiet, comfortable sleep.



For more information, please call us at 1-800-523-7672  
Email: [support@microsonic-inc.com](mailto:support@microsonic-inc.com) or visit: [microsonic-inc.com](http://microsonic-inc.com)  
Woman-owned and operated. ISO 9001:2015 Certified

UNIQUE DESIGN • COMFORTABLE FIT



Patent Pending Design

## LIFEPLUGS®

- **COMFORT** Our super soft silicone is perfect for all day use
- **PROTECTION** NRR 17/25 dB, Mean 35/39 dB
- **HEALTH** Hypoallergenic & Antibacterial
- **STYLE** Thousands of color combinations
- **LIFEPLUGS®** Perfect for yoga, exercise, sports, music concerts, or studying
- **LIFEPLUGS® SLEEP** Perfect for sleeping



## LIFEPLUGS® HEARING PROTECTION SPECIFICATIONS


	LifePlugs®	LifePlugs® Sleep
Use	Hearing Protection	Sleeping
NRR	17	25
Material	Silicone	Silicone
Earplug Color Options	Clear	Clear
	Thousands of Colors	Thousands of Colors
Handle with Lanyard	Cloth or Plastic	
Lanyard Color Options	Black	
	Clear	

## ATTENUATION

ANSI S3.19-1974	Frequency in Hertz								
	125	250	500	1000	2000	3150	4000	6000	8000
Standard	15	14.6	21.9	27.5	34.6	31.6	27.2	27.8	34.2
Sleep	33.3	29.1	33.9	33.2	39.2	37	34.3	31.5	36.8

## DATA

### Attenuation Values for LifePlugs® Custom Hearing Protection Earmolds



TO FIND OUT MORE ABOUT MICROSONIC PLEASE CALL:

800-523-7672 | [support@microsonic-inc.com](mailto:support@microsonic-inc.com) | [microsonic-inc.com](http://microsonic-inc.com)



[www.uvero.com](http://www.uvero.com)

Microsonic Inc. 2960 Duss Ave. Ambridge, PA 15003

Woman-owned & operated. ISO 9001:2015 Certified

©2022 Microsonic Inc. All Rights Reserved.