


Custom earmolds open-fit & RIC/RITE



#OF2	#OF4	#OF5	#OF7	#OF21	#SL2	#SL4	#SL5
OPEN-FIT SKELETON (w/IROS) (w/slim tube)	OPEN-FIT CANAL (w/IROS) (w/slim tube)	OPEN-FIT CANAL-LOK (w/IROS) (w/slim tube)	OPEN-FIT CANAL-SHELL (w/IROS) (w/slim tube)	OPEN-FIT ADV DESIGN FREE FIELD (w/SAV) (w/slim tube)	SLIM-FIT SKELETON Specify vent size & receiver	SLIM-FIT CANAL Specify vent size & receiver	SLIM-FIT CANAL-LOK Specify vent size & receiver
Available in all materials. Provides largest vent possible. SAV is optional	Available in all materials. Provides largest vent possible. SAV is optional	Available in all materials. Provides largest vent possible. SAV is optional	Available in all materials. Provides largest vent possible. SAV is optional	Available in all materials for use with slim tubes. Non-occluding	Recommended for RIC. Soft materials ONLY	Recommended for RIC. Soft materials ONLY	Recommended for RIC. Soft materials ONLY



For more information, please call us at 1-800-523-7672 Email: support@microsonic-inc.com
 Visit: microsonic-inc.com | Woman-owned and operated. ISO 9001:2015 Certified
 ©2022 Microsonic Inc. All Rights Reserved.

