


# MICROSONIC<sup>®</sup>


CLEARLY THE SOUND CHOICE

## Custom hearing solutions HEARING PROTECTION

### PASSGUARD<sup>™</sup>



- NRR 27 dB
- Perfect for high impact noise
- Permits hearing conversation
- Custom made for perfect fit
- Hypoallergenic & Antibacterial
- 3D scanned & archived: Hassle-free re-orders/remakes
- Optional EZ-Grip<sup>™</sup> handles. Cords stay close to ear
- Optional laser engravings: company name, logo, left/right, personal initials, name, serial number, etc.
- Available Styles: Shell, Canal, and Canal-Shell



microsonic-inc.com  
Microsonic Inc.  
CAGE Code: 4DTY2  
2960 Duss Avenue  
Ambridge, PA 15003  
P: 800-523-7672  
F: 724-266-6309  
support@microsonic-inc.com



For more information, please call us at 1-800-523-7672  
Email: [info@microsonic-inc.com](mailto:info@microsonic-inc.com) | Web: [microsonic-inc.com](http://microsonic-inc.com)  
Woman owned and operated. ISO 9001 Certified 2015  
©2024 Microsonic Inc. All Rights Reserved.